

VASCULAR ACCESS

Extended Dwell Peripheral IV Catheter



leaderflex

Reducing venous depletion
for **all** patients



Value Life

When should you use an Extended Dwell Peripheral IV Catheter?

The leaderflex family of 22 Ga catheters are versatile catheters for:

Reduce Repeated Veinpunctures


Vein Preservation

Infusion
(Antibiotics, Hydration & Pain Medications)

Therapy for 1-4 weeks

“ The catheter selected shall be of the **smallest gauge and length** with the **fewest number of lumens** and shall be the **least invasive device** needed to accommodate and manage the prescribed therapy. ”

Infusion Nursing Standards of Practice, 2011:34(15)



• leaderflex

Product code	Catheter (polyurethane)			Puncture needle		Guidewire		Units per box
	Ø (Fr)	Length (cm)	Flow rate (ml/min)	Gauge	Length (cm)	Ø (in)	Length (cm)	
VYLF1004	2	4	17	21	4	0.018	23	10
VYLF1006	2	6	14	21	4	0.018	23	10
VYLF1008	2	8	12	21	7	0.018	26	10
VYLF1020	2	20	4.4	21	7	0.018	50	10

• Accessories

Product code	leaderflex® Accessories and Kits	Units per box
5804.08	Grip-Lok Pediatric Adhesive Securement Device, Individually Packaged	20

For further information, please contact: customerservice@vygonusa.com.

The specifications shown in this brochure are for information only and are not, under any circumstances, of a contractual nature.